

Testing for a COVID-19 Vaccine?

Reserve your Class 2 Biological Safety Cabinet today

AVIDITY
SCIENCE

Avidity Science is prepared to support our key workers in laboratories around the world. Due to the large demand for **Class 2 Biological Safety Cabinets**, we are increasing our stock and you can **reserve yours with us now!**

NuAire Class 2 A2 Recirculating Biological Safety Cabinets

NuAire Class 2 Cabinets meet the BS/EN 12469 and are independently certified by TUV, Technischer Überwachungsverein or Technical Monitoring Association.

- ✓ The German certification board is widely accepted as the most stringent and comprehensive test a safety cabinet can undergo.

NuAire Class 2 Cabinets meet the NSF/ANSI 49 standard

- ✓ The NSF/ANSI 49 standard is regarded as the “gold standard” for biosafety cabinets in the USA.



Why NuAire Cabinets are one of the Safest Class 2 Cabinets available on the market

- ✓ All stainless steel construction, double walled and sealed using welding, producing a one-piece leak proof monolithic structure. Watch this video to learn more: <https://youtu.be/vXPz55pSyhQ>
- ✓ Double exhaust HEPA filters as standard
- ✓ Inflow velocity 0.53m/s provides a safe effective barrier for operator protection
- ✓ Removable single solid work surface makes it easier to clean unlike a segmented surface
- ✓ 3 Litre spill containment and a drain with switch valve
- ✓ A 10 inch access opening and ergonomic armrest provides user comfort
- ✓ Two power sockets on the inside, service Valve RHS wall, service coupling RHS wall and two Service couplings LHS wall, some cabinets come with no service valves
- ✓ Manual frame-less sliding tempered sash window simple and fast with no downtime
- ✓ **Six Years Warranty** includes all parts and ALL HEPA FILTERS that includes the MAIN HEPA and TWO EXHAUST HEPA filters. These are consumables and won't be covered on the other warranty options you may be offered by other companies. If they do, then you should get it in writing as we will on your quote.
- ✓ **Avidity Science** are a Service led company who have 26 Specialist engineers in the UK working locally to support you. We have been installing and servicing NuAire in the UK for over 20 years. Our staff have the experience, dedication and commitment to support you and your vital research.

ORDER AND RESERVE YOURS TODAY - Don't miss out!

For more information or if you require a Quotation, please contact
UK.Info@AvidityScience.com