



## Designed to make working with drinking valves even easier!

The animal drinking valve tools are used to remove a drinking valve from a manifold or cage, disassemble or assemble valves, or adjust the flow setting on the A130 or A140 valve series.



### A160 Valve Tool

P/N: 4500-8781-131  
 TOOL LENGTH: 3.9 in. (99 mm)  
 FITS SERIES: A160  
 BARREL DIAMETER: 0.49 in. (12.4 mm)



### A140 Valve Tool

P/N: 4500-8781-231  
 TOOL LENGTH: 3.2 in. (81 mm)  
 FITS SERIES: Current A140, A141, A145  
 BARREL DIAMETER: 0.23 in. (5.8 mm)



### A142 Valve Tool

P/N: 4500-8781-331  
 TOOL LENGTH: 3.4 in. (86 mm)  
 FITS SERIES: A142, A143  
 BARREL DIAMETER: 0.36 in. (9.1 mm)



### A130 Valve Tool

P/N: 4500-8781-431  
 TOOL LENGTH: 3.2 in. (81 mm)  
 FITS SERIES: Older A130  
 BARREL DIAMETER: 0.25 in. (6.4 mm)



### QD2 3-Slot Tool Bit

HANDLE TYPE: 0.250 in Hex Driver  
 P/N: 4500-8900  
 FITS: QD2 Socket Body of Valves  
 Integral to the Manifold / Series 1660-8900



### RACHETING MULTI-BIT SCREWDRIVER

HANDLE TYPE: 0.250 in Hex Driver  
 P/N: 4500-7011-001  
 FITS: QD2 Socket Body of Valves  
 Integral to the Manifold / Series 1660-8900