



A160 Drinking Valve Identification Chart

A160 Integral to the Manifold



1000-1641
A164 SHIELDED 3/8-32 1.265LG SS



1000-1641-350
A164 SHIELDED VALVE 3/8-32 RAT

Here are some examples of valves with extenders (shown transparent):
Note: Customers to retain the extenders on their manifolds. They should only send in the core '1000-1641' drinking valve for E-vive renewal.

1000-1641



1000-1641-001
A164 SHIELDED W/1.0 EXT LP



1000-1641-002
A164 SHIELDED W/1.5 EXT LP



1000-1641-101
A164 SHIELDED 3.265EXT 02MPS LP RAIR

A160 Quick Disconnect



1000-1661
A166 SHIELDED QD2 1.73 LG



1000-1661-351
A166 SHIELDED QD2 RCSD-STEM 1.80LG



1000-1662
A166 SHIELDED QD2 2.30LG



1000-1662-351
A166 SHIELDED QD2 RCSD-STEM 2.30LG



1000-1665-351
A166 SHIELDED QD2 RCSD-STEM 1.33LG



1000-1665
A166 SHIELDED QD2 W/CV

A160 Quick Disconnect with Cage Seal



1000-1661-001
A166 SHIELDED QD2 1.73"LG W/CAGE-SEAL



1000-1661-350
A166 SHIELDED QD2 RCSD-STEM 1.80IN LG W/SEAL



1000-1662-350
A166 SHIELDED QD2 RCSD-STM 2.30LG W/SEAL



1000-1665-001
A166 SHIELD QD2 W/CAGE SEAL



1000-1665-003
A166 SHIELD QD2 W/VALVE SEAL NATURAL



1000-1665-350
A166 SHIELDED QD2 RCSD-STEM 1.33LG W/SEAL



Note: Customers to retain Tecniplast blue or green cage seals / boots.

See samples in green to the left.



A160 Drinking Valve Identification Chart

A160 Integral to the Cage



1000-1684
A168 SHIELDED EC-LP-SHRT



1000-1688
A168 SHIELDED EC-LP



1000-1689
A168 SHIELDED EC-STND



1000-1684-350
A168 SHIELDED EC-LP-SHRT RCSD-STEM




1000-1688-350
A168 SHIELDED EC-LP RCSD-STEM



1000-1689-350
A168 SHIELDED EC-STND RCSD-STEM

Comments / Notes

- Customer's with A160 valve assemblies that are integral to the cage should be advised to retain the accompanying hardware (IE - nut, washers & seals) and send only the valve for E-vive care. See diagram below.



Sample valve retaining hardware used in assemblies integral to the cage.

Valve Part Number

1665
03312009

Manufacture Date

Silicone Shield Color also Indicates Date of Manufacture



GREEN	RED	GRAY	ORANGE	BROWN	YELLOW	BLUE	PINK	TEAL
Only Used From 2000 - 2006	Only Used In 2007	2008, 2015	2009, 2016	2010, 2017	2011, 2018	2012, 2019	2013, 2020	2014, 2021